



Insight into health

## 3000K

Full Digital Ultrasound System

---



## 3000K

Full Digital Ultrasound System

---

**Wuxi Biomedical Technology Co., Ltd.**

Add: Room 408, No. 801 Hongqiao road, Wuxi, Jiangsu,  
China 214072  
Tel: 86-510-85572789  
Fax: 86-510-85572789  
Website: <http://www.portable-ultrasound-scanner.com>

# 3000K

## Full Digital Ultrasound System

**Application**  
 Abdomen OB&GYN Urology Blood vessel Small Organ  
 Pediatric Musculoskeletal Basic Cardiac

BIO 3000K is a notebook Black and White ultrasound system with full touch screen panel which reflects BIO ultimate pursuit of ultrasound image, operation and mobile performance. With a break through design concept, we provide a more clear image quality, more stable performance, and bring you a comprehensive clinical diagnosis solution, giving you a new experience for your ultrasonic diagnostic work.



- 9.7 inch touch screen control panel makes the operation of the system more simple
- Built-in 32GB SSD for more reliable data storage
- Built-in lithium battery for 2 hours of continuous operation
- 12 inch LED display

### Optional Transducer

Probe Model	Probe type	Frequency Range
CR60	Convex probe	2-5Mhz
MCR11	Micro-convex probe	5-8Mhz
MCR20	Micro-convex probe	2-4Mhz
L25	Linear probe	6-15Mhz
L40	Linear probe	5-12Mhz
MCVR10	Transvaginal probe	5-8Mhz

### Advanced Technology

- **Efficiency**  
 High efficient workflow design, allows the doctors complete the daily routine inspection quickly and easily, which greatly reduces their workload.
- **Intelligent**  
 Intelligent image optimization solutions can help you carry out diagnosis and assessment of the difficult patients with more confidence. It makes your operation more accurate and efficient.
- **Visualization**  
 Portable notebook design, full touch screen control panel makes machine interaction more intuitive, convenient, and allows doctors to master visual patient information anytime, anywhere.

